

Physician's Prescription for Mechanical Ventilation with the LTV 1150

Patient Name:			
Patient DOB:			
Front Panel Settings Physician Must Order			
Breath type:	<input type="checkbox"/> Volume <u>or</u> <input type="checkbox"/> Pressure	Mode:	<input type="checkbox"/> Assist/Control <input type="checkbox"/> SIMV / CPAP <input type="checkbox"/> NPPV*
Breath Rate:		Tidal Volume:	
Pressure Control*:		Inspiratory Time:	
Pressure Support*:		O ₂ % (lpm):	
PEEP:		*NPPV Settings:	<input type="checkbox"/> IPAP <input type="checkbox"/> EPAP

**LTV 1150 is PEEP Compensated. So, if set pressure at 10 cm, and set PEEP of 5 cm, the PIP will be 15 cm.*

The remaining settings may be ordered by the physician, or the physician may elect for the RT to determine the appropriate settings.			
Flow Trigger Sensitivity Setting – Front Panel			
<input type="checkbox"/> RT to determine sensitivity <u>or</u> <input type="checkbox"/> Physician orders: _____			
Alarm Settings – Front Panel			
<input type="checkbox"/> RT to set the alarm-related front panel settings <u>or</u> <input type="checkbox"/> Physician orders alarm settings as indicated below:			
High Pressure Limit:	<input type="checkbox"/> PIP + 10 cm <input type="checkbox"/>	Low Pressure Limit:	<input type="checkbox"/> PIP – 5 cm (not < 5 cm) <input type="checkbox"/>
Low Minute Volume:	<input type="checkbox"/> approximately 25% less than patient's sleeping VE <input type="checkbox"/>		
Alarm Operations – Extended Menu			
<input type="checkbox"/> RT to set the alarm-related extended menu items <u>or</u> <input type="checkbox"/> Physician orders alarm-related extended menu items as indicated below:			
Alarm Volume:	<input type="checkbox"/> 85 decibels <input type="checkbox"/>	Apnea Interval:	<input type="checkbox"/> 20 seconds <input type="checkbox"/>
Hi Pressure Delay:	<input type="checkbox"/> No delay <input type="checkbox"/> 1 breath <input type="checkbox"/> 2 breaths	LPP Alarm:	<input type="checkbox"/> All breaths <input type="checkbox"/> PC/VC breaths only
Hi PEEP:	<input type="checkbox"/> 5 above set PEEP <input type="checkbox"/>	Hi Rate & Time before alarm:	<input type="checkbox"/> Preset setting (40, 60 or 80) <input type="checkbox"/> <input type="checkbox"/> Time: 30 sec (preset setting) <input type="checkbox"/> Time:
Low PEEP:	<input type="checkbox"/> -3 below set PEEP <input type="checkbox"/>		
Vent Operations – Extended Menu			
<input type="checkbox"/> RT to set the vent operations extended menu items <u>or</u> <input type="checkbox"/> Physician orders vent operations extended menu items as indicated below:			
Rise Time:	<input type="checkbox"/>	Flow Term:	<input type="checkbox"/>
Time Term:	<input type="checkbox"/> 1.5 X IT <input type="checkbox"/>	PC Flow Term:	<input type="checkbox"/> On <input type="checkbox"/> Off
Control Unlock:	<input type="checkbox"/> Easy <input type="checkbox"/> Hard		

Physician Signature _____ Date _____